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Article Clinisafe 01315 - Nitrile powder free glove, single use, non sterile

Informative Note of Protection Glove against chemical risks and microorganisms. Caution: Read carefully before use.

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Pre-use Checks and Instruction for Use.

Before the use, check visually the PPE to be sure of its integrity. If a glove is damaged, it should not be used. This PPE is disposable. Do not re-use it.

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Labeling

The CE mark shows that the glove is in conformity with the essential requirements of Regulation EU 2016/425 on Personal Protective Equipment, such as ergonomics, safe, comfort, and that it is produced with materials that are not based on or are not treated with known allergic substances. This device is certified according to EN 420:2003 +A1 2009, EN ISO 374-1:2016, EN 374-2:2014; EN 374-4:2013 e EN ISO 374-5:2016.

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Protection against chemical agents

This personal protective equipment PPE is designed to protect the user against accidental contact with some chemicals and micro-organism that could cause damage to the user himself in the exercise of his work.

The level of performance is determined by the result of the correspondent test as shown in EN ISO 374-1:2016 standard for protection against chemicals.

The levels of performance are based on the test results done in laboratory: thus, they are not necessarily equal to the work environmental existing conditions. The protection level gives the time necessary to the permeation of the chemical tested. The levels of performance indicated as per EN ISO 374-1:2016 standard, show the ability of the glove to protect against one specific agent for a defined period of time of potential accidental contact with, as described in the following table:

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This PPE offers the following levels of protection, determined on the laboratory test results in conformity with EN 16523-1:2015 standard for the resistance to permeation and in conformity with EN 374-4: 2013 standard for degradation resistance as indicated in EN ISO 374-1:2016 standard. The test chemicals and the correspondent levels of protection are shown in table A.

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This PPE is classified as Type B because it has a level 2 permeation level for 3 chemical compounds (G-K -J) as per table 2 EN ISO 374-1:2016 standard.

Is not offered any protection against any other kind of potential risk that is non-explicitly indicated.

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Protection against biological risks

The glove was submitted to the penetration test in conformity with EN 374-2:2014 standard, and it passed the water leak test and the air leak test. It is conventionally assumed that, offering a good barrier against water and air leaks, the PPE offers protection against the biological agents: bacteria and fungi. This evaluation is, however, a conventional conclusion.

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Storage

The glove should be stored avoiding the direct lights and the solar radiations, at a environment temperature between 5°C and 40°C. Store in a cool, dry place. Avoid ozone. Remove the PPE from the box just before the use.

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Cleaning

Not applicable. The PPE is disposable.

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Disposal

The product must be disposed of in compliance with the applicable local regulations. If the glove has been used in contact with biological fluids, they must be disposed of in the appropriate containers for hazardous and infectious waste.

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Validity of the product

The validity period of the product is five years from the production date.

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This PPE was certified by: A.N.C.I. Servizi Sezione C.I.M.A.C. Corso Brodolini, 19 - 27029 Vigevano (PV) – Italy – identification number: 0465.

The Responsible of the Prevention and Protection Department shall inform the user on the content of this informative note.